

# MOGT3 Polyclonal Antibody

Cat No: HR1AP10405

For research use only

## Overview

Product Name	MOGT3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	MOGAT3 DC7 DGAT2L7 UNQ9383/PRO34208
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	monoacylglycerol O-acyltransferase 3(MOGAT3) Homo sapiens Acyl-CoA:monoacylglycerol acyltransferase (MOGAT; EC 2.3.1.22) catalyzes the synthesis of diacylglycerol from 2-monoacylglycerol and fatty acyl-CoA (Cheng et al., 2003 [PubMed 12618427]).[supplied by OMIM, Mar 2008],