

## PO5F2 Polyclonal Antibody

Cat No: HR1AP10396

For research use only

### Overview

Product Name	PO5F2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	POU5F2 SPRM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Transcription factor that binds preferentially to the octamer motif (5'-ATGTTAAT-3'). May exert a regulatory function in meiotic events that are required for terminal differentiation of male germ cell.,similarity:Belongs to the POU transcription factor family. Class-5 subfamily.,similarity:Contains 1 homeobox DNA-binding domain.,similarity:Contains 1 POU-specific domain.,tissue specificity:Expressed in skeletal and cardiac muscles, brain, heart and lung. Little or no detectable expression found in pancreas, kidney, liver or placenta.,</p>