

PO3F3 Polyclonal Antibody

Cat No: HR1AP10395

For research use only

Overview

Product Name	PO3F3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	POU3F3 BRN1 OTF8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	POU class 3 homeobox 3(POU3F3) Homo sapiens This gene encodes a POU-domain containing protein that functions as a transcription factor. The encoded protein recognizes an octamer sequence in the DNA of target genes. This protein may play a role in development of the nervous system. [provided by RefSeq, Apr 2015],