

FMR1N Polyclonal Antibody

Cat No: HR1AP10391

For research use only

Overview

Product Name	FMR1N Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	FMR1NB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 P-type (trefoil) domain.,tissue specificity:Testis-specific. Expressed in melanoma, sarcoma, lung, breast, bladder, esophageal and ovarian cancers.,