

TNIP3 Polyclonal Antibody

Cat No: HR1AP10390

For research use only

Overview

Product Name	TNIP3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	TNIP3 ABIN3 LIND
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	induction:By Listeria infection. Expression is slightly down-regulated by dexamethasone and slightly up-regulated by IL-10.,tissue specificity:Highly expressed in lung, lymph node, thymus and fetal liver. Expressed at lower levels in brain, spleen, leukocytes and tonsils.,