

SATB2 Polyclonal Antibody

Cat No: HR1AP10386

For research use only

Overview

Product Name	SATB2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	80kDa
GeneID?Human?	SATB2 KIAA1034
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	SATB homeobox 2(SATB2) Homo sapiens This gene encodes a DNA binding protein that specifically binds nuclear matrix attachment regions. The encoded protein is involved in transcription regulation and chromatin remodeling. Defects in this gene are associated with isolated cleft palate and mental retardation. Alternate splicing results in multiple transcript variants that encode the same protein. [provided by RefSeq, Feb 2010].