

## BCL6B Polyclonal Antibody

Cat No: HR1AP10381

For research use only

### Overview

Product Name	BCL6B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	BCL6B BAZF ZNF62
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Acts as a sequence-specific transcriptional repressor in association with BCL6. May function in a narrow stage or be related to some events in the early B-cell development.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 C2H2-type zinc fingers.,subunit:Associates with BCL6 through the BTB domain.,tissue specificity:Ubiquitously expressed with higher expression found in heart and placenta.,