

ALKB2 Polyclonal Antibody

Cat No: HR1AP10379

For research use only

Overview

Product Name	ALKB2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	ALKBH2 ABH2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>alkB homolog 2, alpha-ketoglutarate dependent dioxygenase(ALKBH2) Homo sapiens The Escherichia coli AlkB protein protects against the cytotoxicity of methylating agents by repair of the specific DNA lesions generated in single-stranded DNA. ALKBH2 and ALKBH3 (MIM 610603) are E. coli AlkB homologs that catalyze the removal of 1-methyladenine and 3-methylcytosine (Duncan et al., 2002 [PubMed 12486230]).[supplied by OMIM, Mar 2008],</p>