

GOT1B Polyclonal Antibody

Cat No: HR1AP10377

For research use only

Overview

Product Name	GOT1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	15kDa
GeneID?Human?	GOLT1B GOT1A CGI-141 HDCMA39P UNQ432/PRO793
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1 family.,