

CR3L4 Polyclonal Antibody

Cat No: HR1AP10376

For research use only

Overview

Product Name	CR3L4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	CREB3L4 AIBZIP CREB4 JAL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cAMP responsive element binding protein 3 like 4(CREB3L4) Homo sapiens This gene encodes a CREB (cAMP responsive element binding) protein with a transmembrane domain which localizes it to the ER membrane. The encoded protein is a transcriptional activator which contains a dimerization domain, and this protein may function in a number of processing pathways including protein processing. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011],</p>