

IKBL1 Polyclonal Antibody

Cat No: HR1AP10372

For research use only

Overview

Product Name	IKBL1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	41kDa
GeneID?Human?	NFKBIL1 IKBL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NFKB inhibitor like 1(NFKBIL1) Homo sapiens This gene encodes a divergent member of the I-kappa-B family of proteins. Its function has not been determined. The gene lies within the major histocompatibility complex (MHC) class I region on chromosome 6. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2009],