

# STRN4 Polyclonal Antibody

Cat No: HR1AP10371

For research use only

## Overview

Product Name	STRN4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	82kDa
GeneID?Human?	STRN4 ZIN
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function: Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein., miscellaneous: The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane)., similarity: Belongs to the WD repeat striatin family., similarity: Contains 7 WD repeats.,