

THIM Polyclonal Antibody

Cat No: HR1AP10367

For research use only

Overview

Product Name	THIM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	ACAA2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	acetyl-CoA acyltransferase 2(ACAA2) Homo sapiens The encoded protein catalyzes the last step of the mitochondrial fatty acid beta-oxidation spiral. Unlike most mitochondrial matrix proteins, it contains a non-cleavable amino-terminal targeting signal. [provided by RefSeq, Jul 2008],