

CENPK Polyclonal Antibody

Cat No: HR1AP10364

For research use only

Overview

Product Name	CENPK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	CENPK ICEN37 FKSG14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	centromere protein K(CENPK) Homo sapiens CENPK is a subunit of a CENPH (MIM 605607)-CENPI (MIM 300065)-associated centromeric complex that targets CENPA (MIM 117139) to centromeres and is required for proper kinetochore function and mitotic progression (Okada et al., 2006 [PubMed 16622420]).[supplied by OMIM, Mar 2008],