

CENPC Polyclonal Antibody

Cat No: HR1AP10362

For research use only

Overview

Product Name	CENPC Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	103kDa
GeneID?Human?	CENPC1 CENPC ICEN7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	centromere protein C(CENPC) Homo sapiens Centromere protein C 1 is a centromere autoantigen and a component of the inner kinetochore plate. The protein is required for maintaining proper kinetochore size and a timely transition to anaphase. A putative pseudogene exists on chromosome 12. [provided by RefSeq, Jul 2008],