

TNIK Polyclonal Antibody

Cat No: HR1AP10356

For research use only

Overview

Product Name	TNIK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	149kDa
GeneID?Human?	TNIK KIAA0551
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TRAF2 and NCK interacting kinase(TNIK) Homo sapiens Germinal center kinases (GCKs), such as TNIK, are characterized by an N-terminal kinase domain and a C-terminal GCK domain that serves a regulatory function (Fu et al., 1999 [PubMed 10521462]).[supplied by OMIM, Mar 2008],