

CRADD Polyclonal Antibody

Cat No: HR1AP10353

For research use only

Overview

Product Name	CRADD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	CRADD RAIDD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>CASP2 and RIPK1 domain containing adaptor with death domain(CRADD) Homo sapiens This gene encodes a protein containing a death domain (DD) motif. This protein recruits caspase 2/ICH1 to the cell death signal transduction complex, which includes tumor necrosis factor receptor 1 (TNFR1A) and RIPK1/RIP kinase, and acts in promoting apoptosis. A mutation in this gene was associated with mental retardation. A related pseudogene is found on chromosome 3. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016],</p>