

# MLRS Polyclonal Antibody

Cat No: HR1AP10351

For research use only

## Overview

Product Name	MLRS Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	MYLPF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	miscellaneous:This chain binds calcium.,similarity:Contains 3 EF-hand domains.,subunit:Myosin is a hexamer of 2 heavy chains and 4 light chains.,tissue specificity:Expressed in fetal and adult skeletal muscle.,