

WNT7B Polyclonal Antibody

Cat No: HR1AP10345

For research use only

Overview

Product Name	WNT7B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	WNT7B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>Wnt family member 7B(WNT7B) Homo sapiens This gene is a member of the WNT gene family, which consists of structurally related genes that encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. Among members of the human WNT family, this gene product is most similar to WNT7A protein. [provided by RefSeq, Oct 2008],</p>