

VEZF1 Polyclonal Antibody

Cat No: HR1AP10338

For research use only

Overview

Product Name	VEZF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	VEZF1 DB1 ZNF161
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	vascular endothelial zinc finger 1(VEZF1) Homo sapiens Transcriptional regulatory proteins containing tandemly repeated zinc finger domains are thought to be involved in both normal and abnormal cellular proliferation and differentiation. ZNF161 is a C2H2-type zinc finger protein (Koyano-Nakagawa et al., 1994 [PubMed 8035792]). See MIM 603971 for general information on zinc finger proteins.[supplied by OMIM, Sep 2008],