

TIAF1 Polyclonal Antibody

Cat No: HR1AP10337

For research use only

Overview

Product Name	TIAF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	TIAF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:The TIAF1 protein is coded in the 3'-UTR region of MYO18A.,function:Inhibits the cytotoxic effects of TNF-alpha and overexpressed TNF receptor adapters TRADD, FADD, and RIP. Involved in TGF-beta1 inhibition of IkappaB-alpha expression and suppression of TNF-mediated IkappaB-alpha degradation.,