

## PLCH2 Polyclonal Antibody

Cat No: HR1AP10334

For research use only

### Overview

Product Name	PLCH2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	155kDa
GeneID?Human?	PLCH2 KIAA0450 PLCL4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	phospholipase C eta 2(PLCH2) Homo sapiens PLCH2 is a member of the PLC-eta family of the phosphoinositide-specific phospholipase C (PLC) superfamily of enzymes that cleave PtdIns(4,5) P2 to generate second messengers inositol 1,4,5-trisphosphate and diacylglycerol (Zhou et al., 2005 [PubMed 16107206]).[supplied by OMIM, Jun 2009],