

PPBN Polyclonal Antibody

Cat No: HR1AP10335

For research use only

Overview

Product Name	PPBN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GeneID?Human?	ALPPL2 ALPPL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	alkaline phosphatase, placental like 2(ALPPL2) Homo sapiens There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The product of this gene is a membrane bound glycosylated enzyme, localized to testis, thymus and certain germ cell tumors, that is closely related to both the placental and intestinal forms of alkaline phosphatase. [provided by RefSeq, Jul 2008],