

PDD2L Polyclonal Antibody

Cat No: HR1AP10332

For research use only

Overview

Product Name	PDD2L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	PDCD2L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	PDCD2L (Programmed Cell Death 2 Like) is a Protein Coding gene. An important paralog of this gene is ENSG00000266953. Over-expression suppresses AP1, CREB, NFAT, and NF-κB transcriptional activation, and delays cell cycle progression at S phase.