

## PDCD7 Polyclonal Antibody

Cat No: HR1AP10331

For research use only

### Overview

Product Name	PDCD7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	PDCD7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	programmed cell death 7(PDCD7) Homo sapiens This gene encodes a 59 kDa protein that is associated with the U11 small nuclear ribonucleoprotein (snRNP), which is a component of the minor U12-type spliceosome responsible for catalyzing pre-mRNA splicing of U12-type introns. [provided by RefSeq, Dec 2010],