

NFAM1 Polyclonal Antibody

Cat No: HR1AP10327

For research use only

Overview

Product Name	NFAM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	NFAM1 CNAIP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	NFAT activating protein with ITAM motif 1(NFAM1) Homo sapiens The protein encoded by this gene is a type I membrane receptor that activates cytokine gene promoters such as the IL-13 and TNF-alpha promoters. The encoded protein contains an immunoreceptor tyrosine-based activation motif (ITAM) and is thought to regulate the signaling and development of B-cells. [provided by RefSeq, Jul 2008],