

## NB5R4 Polyclonal Antibody

Cat No: HR1AP10326

For research use only

### Overview

Product Name	NB5R4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	CYB5R4 NCB5OR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cytochrome b5 reductase 4(CYB5R4) Homo sapiens NCB5OR is a flavohemoprotein that contains functional domains found in both cytochrome b5 (CYB5A; MIM 613218) and CYB5 reductase (CYB5R3; MIM 613213) (Zhu et al., 1999 [PubMed 10611283]).[supplied by OMIM, Jan 2010],