

# MYL6B Polyclonal Antibody

Cat No: HR1AP10323

For research use only

## Overview

Product Name	MYL6B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	MYL6B MLC1SA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>myosin light chain 6B(MYL6B) Homo sapiens Myosin is a hexameric ATPase cellular motor protein. It is composed of two heavy chains, two nonphosphorylatable alkali light chains, and two phosphorylatable regulatory light chains. This gene encodes a myosin alkali light chain expressed in both slow-twitch skeletal muscle and in nonmuscle tissue. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2010],</p>