

MYLK2 Polyclonal Antibody

Cat No: HR1AP10324

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | MYLK2 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Mouse,Rat |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 65kDa |
| GeneID?Human? | MYLK2 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | myosin light chain kinase 2(MYLK2) Homo sapiens This gene encodes a myosin light chain kinase, a calcium/calmodulin dependent enzyme, that is exclusively expressed in adult skeletal muscle. [provided by RefSeq, Jul 2008], |