

MAGF1 Polyclonal Antibody

Cat No: HR1AP10318

For research use only

Overview

Product Name	MAGF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	MAGEF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MAGE family member F1(MAGEF1) Homo sapiens This intronless gene encodes a member of the MAGE superfamily. It is ubiquitously expressed in normal tissues and in tumor cells. This gene includes a microsatellite repeat in the coding region. [provided by RefSeq, Jul 2008],