

MAGC3 Polyclonal Antibody

Cat No: HR1AP10316

For research use only

Overview

Product Name	MAGC3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	MAGEC3 HCA2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MAGE family member C3(MAGEC3) Homo sapiens This gene is a member of the MAGEC gene family. The members of this family are not expressed in normal tissues, except for testis, and are expressed in tumors of various histological types. The MAGEC genes are clustered on chromosome Xq26-q27. Two transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008],