

MAGB5 Polyclonal Antibody

Cat No: HR1AP10315

For research use only

Overview

Product Name	MAGB5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	MAGEB5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 MAGE domain.,tissue specificity:Expressed in testis. Not expressed in other normal tissues, but is expressed in tumors of different histological origins.,