

K1H1 Polyclonal Antibody

Cat No: HR1AP10301

For research use only

Overview

Product Name	K1H1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	KRT31 HHA1 HKA1 KRTHA1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	keratin 31(KRT31) Homo sapiens The protein encoded by this gene is a member of the keratin gene family. As a type I hair keratin, it is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. [provided by RefSeq, Jul 2008],