



K1C9 Polyclonal Antibody

Cat No: HR1AP10300

For research use only

Overview

Product Name	K1C9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	68kDa
GeneID?Human?	KRT9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	keratin 9(KRT9) Homo sapiens This gene encodes the type I keratin 9, an intermediate filament chain expressed only in the terminally differentiated epidermis of palms and soles. Mutations in this gene cause epidermolytic palmoplantar keratoderma. [provided by RefSeq, Jul 2008],