

# K1C12 Polyclonal Antibody

Cat No: HR1AP10298

For research use only

## Overview

Product Name	K1C12 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	KRT12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	keratin 12(KRT12) Homo sapiens KRT12 encodes the type I intermediate filament chain keratin 12, expressed in corneal epithelia. Mutations in this gene lead to Meesmann corneal dystrophy. [provided by RefSeq, Jul 2008],