

ITA11 Polyclonal Antibody

Cat No: HR1AP10292

For research use only

Overview

Product Name	ITA11 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	130kDa
GeneID?Human?	ITGA11 MSTP018
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>integrin subunit alpha 11(ITGA11) Homo sapiens This gene encodes an alpha integrin. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This protein contains an I domain, is expressed in muscle tissue, dimerizes with beta 1 integrin in vitro, and appears to bind collagen in this form. Therefore, the protein may be involved in attaching muscle tissue to the extracellular matrix. Alternative transcriptional splice variants have been found for this gene, but their biological validity is not determined. [provided by RefSeq, Jul 2008].</p>