



INSL5 Polyclonal Antibody

Cat No: HR1AP10291

For research use only

Overview

Product Name	INSL5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	INSL5 UNQ156/PRO182
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	insulin like 5(INSL5) Homo sapiens The protein encoded by this gene contains a classical signature of the insulin superfamily and is highly similar to relaxin 3 (RLN3/INSL7). [provided by RefSeq, Jul 2008],