

INGR2 Polyclonal Antibody

Cat No: HR1AP10288

For research use only

Overview

Product Name	INGR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	IFNGR2 IFNGT1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>interferon gamma receptor 2 (interferon gamma transducer 1)(IFNGR2) Homo sapiens This gene (IFNGR2) encodes the non-ligand-binding beta chain of the gamma interferon receptor. Human interferon-gamma receptor is a heterodimer of IFNGR1 and IFNGR2. Defects in IFNGR2 are a cause of mendelian susceptibility to mycobacterial disease (MSMD), also known as familial disseminated atypical mycobacterial infection. MSMD is a genetically heterogeneous disease with autosomal recessive, autosomal dominant or X-linked inheritance. [provided by RefSeq, Jul 2008].</p>