

# IFI27 Polyclonal Antibody

Cat No: HR1AP10289

For research use only

## Overview

Product Name	IFI27 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	13kDa
GeneID?Human?	IFI27
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	induction:By alpha interferons.,similarity:Belongs to the IFI6/IFI27 family.,