

## IIGP5 Polyclonal Antibody

Cat No: HR1AP10285

For research use only

### Overview

Product Name	IIGP5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	IRGC IIGP5 IRGC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	catalytic activity:GTP + H(2)O = GDP + phosphate.,similarity:Belongs to the interferon-inducible GTPase family.,