

IFT1B Polyclonal Antibody

Cat No: HR1AP10283

For research use only

Overview

Product Name	IFT1B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	IFIT1B IFIT1L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,