

IFRD2 Polyclonal Antibody

Cat No: HR1AP10282

For research use only

Overview

Product Name	IFRD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	IFRD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the IFRD family.,tissue specificity:Expressed in a variety of tissues.,