

IFNK Polyclonal Antibody

Cat No: HR1AP10281

For research use only

Overview

Product Name	IFNK Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	IFNK UNQ6124/PRO20084
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	interferon kappa(IFNK) Homo sapiens This gene encodes a member of the type I interferon family. Type I interferons are a group of related glycoproteins that play an important role in host defenses against viral infections. This protein is expressed in keratinocytes and the gene is found on chromosome 9, adjacent to the type I interferon cluster. [provided by RefSeq, Jul 2008],