

IFN21 Polyclonal Antibody

Cat No: HR1AP10276

For research use only

Overview

Product Name	IFN21 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	IFNA21
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>interferon alpha 21 (IFNA21) Homo sapiens This gene is a member of the alpha interferon gene cluster on the short arm of chromosome 9. Interferons are cytokines produced in response to viral infection that mediate the immune response and interfere with viral replication. The encoded protein is a type I interferon and may play a specific role in the antiviral response to rubella virus. [provided by RefSeq, Sep 2011],</p>