

IFN10 Polyclonal Antibody

Cat No: HR1AP10272

For research use only

Overview

Product Name	IFN10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	IFNA10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>interferon alpha 10(IFNA10) Homo sapiens This gene encodes a protein that belongs to the type I interferon family of proteins, and is located in a cluster of alpha interferon genes on chromosome 9. Interferons are small regulatory molecules that function in cell signaling in response to viruses and other pathogens or tumor cells. This gene is intronless and the encoded protein is secreted. [provided by RefSeq, Aug 2013],</p>