

IFIT3 Polyclonal Antibody

Cat No: HR1AP10270

For research use only

Overview

Product Name	IFIT3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	IFIT3 CIG-49 IFI60 IFIT4 ISG60
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,