

IFIT2 Polyclonal Antibody

Cat No: HR1AP10269

For research use only

Overview

Product Name	IFIT2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	IFIT2 CIG-42 G10P2 IFI54 ISG54
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	induction:By interferons.,similarity:Belongs to the IFIT family.,similarity:Contains 6 TPR repeats.,