

IFI44 Polyclonal Antibody

Cat No: HR1AP10268

For research use only

Overview

Product Name	IFI44 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	IFI44 MTAP44
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:This protein aggregates to form microtubular structures.,induction:By alpha and beta interferons, but not by gamma interferons.,similarity:Belongs to the IFI44 family.,