

IF2B1 Polyclonal Antibody

Cat No: HR1AP10267

For research use only

Overview

Product Name	IF2B1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GeneID?Human?	IGF2BP1 CRDBP VICKZ1 ZBP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>insulin like growth factor 2 mRNA binding protein 1(IGF2BP1) Homo sapiens This gene encodes a member of the insulin-like growth factor 2 mRNA-binding protein family. The protein encoded by this gene contains four K homology domains and two RNA recognition motifs. It functions by binding to the mRNAs of certain genes, including insulin-like growth factor 2, beta-actin and beta-transducin repeat-containing protein, and regulating their translation. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2009],</p>