

IBP5 Polyclonal Antibody

Cat No: HR1AP10266

For research use only

Overview

Product Name	IBP5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	IGFBP5 IBP5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.,similarity:Contains 1 IGFBP N-terminal domain.,similarity:Contains 1 thyroglobulin type-1 domain.,tissue specificity:Osteosarcoma, and at lower levels in liver, kidney and brain.,